INSTRUCTIONS FOR USE – PRECORDIOR AF APP

General

The Precordior AF App mobile application is intended for detection of atrial fibrillation. The measurement follows the activity of the heart based upon movements of the chest.

Measurement

Please remove any thick clothing from your chest. Lie down on your back. Take a comfortable position, as the measurement takes about one minute. Start the measurement by pressing the Start button, which is followed by a 10-second countdown. During the countdown, place the device on your chest. During the measurement, beeps indicate the progress of the measurement event. You will know that the measurement is complete when you hear an audible alarm. You can then remove the phone from your chest. Interpretation of measurement data starts immediately upon completion of the measurement. Depending on the speed of your network connection, results will be available within approximately one minute after the end of the measurement event.

Results

The measurement result indicates your heart rate and whether or not signs of atrial fibrillation were detected. If there is a problem during the measurement, it may be interrupted prematurely. In this case, you will be requested to take a new measurement. In the case of a result indicating signs of atrial fibrillation, you will also be asked to take a second measurement. Results will not be automatically conveyed to third parties. Therefore it is important that the user takes responsibility for possible follow-up actions with healthcare professionals.

Measurement history

All measurement results are automatically stored in the measurement history. If necessary, individual results can be removed from the measurement history view.

Further information on our service

Do you have any questions about our service? Our customer service is available at info@precordior.com.

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Publisher: Precordior Oy, Finland.
Julkaisija: Precordior Oy, Suomi.